

# Reagent Preparation Table

for the *CLA<sup>®</sup> Allergy Test*

The following tables are intended for use with variable-volume pipettes. For greatest accuracy, pipettes with narrow volume ranges should be used.

<b>CLA<sup>®</sup> PETTES</b>	<b>PHOTOREAGENTS</b>
Number per Run	Amount of Each (mL)
1	0.4
2	0.8
3	1.2
4	1.5
5	2.0
6	2.5
7	3.0
8	3.0
9	3.5
10	4.0
11	4.5
12	4.5
13	5.0
14	5.5
15	6.0
16	6.0
17	6.5
18	7.0
19	7.5
20	8.0

In general, photoreagent volumes are rounded to the nearest 0.5 mL, and provide total volumes of at least 1.5 mL per *Pette*. Volumes greater than 1.0 mL may be measured with 1.0 - 5.0 mL variable-volume pipettes, or with 5.0 or 10.0 mL disposable pipettes. If reusable pipettes are used they must be rinsed thoroughly with distilled or de-ionized water to assure that no traces of chemicals, detergents or other contaminants remain.

# HITACHI

#### European Office

Hitachi Chemical Diagnostics, Inc.  
Hitachi Europe, Ltd.  
Whitebrook Park  
Lower Cookham Road  
Maidenhead, Berkshire SL6 8YA  
United Kingdom  
Tel: 44 (0) 1628 585 590

#### United States Office

Hitachi Chemical Diagnostics, Inc.  
630 Clyde Court  
Mountain View, CA 94043-2239  
Tel: 800 233 6278  
[www.hcdiagnostics.com](http://www.hcdiagnostics.com)

Doc. No. 0528 Rev. 06