HITACHI

 Hitachi Chemical Diagnostics, Inc.

CLA[®] IgE ALLERGY POSITIVE CONTROL REAGENT

Doc. No. 0442 - ENG
Rev.: 18
Lot No.: 38198652
Control Kit Part Number: 82060
Bottle Part Number: 93522
Control Kit Contents: 3 x 3 mL

1 Intended Use

CLA IgE Positive Control Reagent is intended for use in quality control procedures, which evaluate performance of the CLA Allergen-Specific IgE Assay.

2 Summary and Principle

The use of quality control materials is indicated as an objective assessment of assay performance. This positive control has been assayed for various specific allergens and their expected values are provided. Allergens that are not listed here may test positive with this lot of control reagent.

3 Reagent

CLA IgE Positive Control Reagent is prepared from human serum. The control is provided in liquid form.

4 Storage and Stability

When stored at $-20 \pm 10^{\circ}\text{C}$, the CLA IgE Positive Control Reagent is stable until the expiration date indicated on the vial label. The reagent is supplied in volumes of 3 mL. When ready to use, thaw, remove the volume needed and re-freeze the unused material immediately. Do not store in the refrigerator. Unused portions of control reagents may be frozen at $-20 \pm 10^{\circ}\text{C}$. Only one freeze-thaw cycle is recommended. Thus aliquoting the control reagents into smaller volumes before refreezing may be helpful.

5 Expected Values

The expected ranges printed in this insert were established by Hitachi Chemical Diagnostics, Inc. and are specific to this lot of CLA IgE Positive Control Reagent. Individual laboratory means should fall within the expected range for each allergen. Procedural or reagent modifications may alter the mean values. Each laboratory should establish its own mean values and acceptable ranges and use those provided only as guides. Hitachi Chemical Diagnostics, Inc. defines acceptable performance criteria for this control as 85% of the assayed target allergens falling within the expected class range.

6 Procedure

The CLA IgE Positive Control Reagent should be tested using the same procedure as patient specimens. Frequency of use may be determined by each Laboratory's Quality Control policies. Allow control reagent to equilibrate to room temperature. Invert gently several times and centrifuge for 10-20 minutes at 2000-3000xg or 2500-3600 rpm immediately prior to use.

7 Limitations

The CLA IgE Positive Control Reagent should not be used past the expiration date. If there is evidence of microbial contamination or excessive turbidity, discard the vial.

WARNING

This product contains human source material. Treat as potentially infectious.


Each donor unit used to manufacture this control was tested by FDA accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Antigen, and antibody to HIV-1. No test method can offer complete assurance that products containing human source materials will be absent of these and other infectious agents. Handle this product with the same precautions used with patient specimens.

FOR IN VITRO DIAGNOSTIC USE ONLY

Allergens	Expected Values	
	Class Range	NCCLS Code
Alternaria	1-3	m6
Aspergillus	3-4	m3
Avocado	1-3	f96
Barley	1-3	f6
Black Alder	2-4	t2
Carrot	1-3	f31
Cat	3-4	e1
Cladosporium	1-3	m2
Cultivated Rye	3-4	g12
Dog	2-4	e5
Egg White	2-4	f1
English Plantain	1-3	w9
Garlic	2-4	f47
Horse	1/0-2	e3
Kentucky Blue Grass	3-4	g8
London Planetree	1-3	t11
Mite, farinae	2-4	d2
Mite, pteronyssinus	2-4	d1
Mugwort	1-3	w6
Oat	2-4	g14
Olive Tree	1/0-2	t9
Orchard Grass	3-4	g3
Pellitory	1-3	w21
Penicillium	2-4	m1
Pork	1/0-2	f26
Rye	1-3	f5
Short Ragweed	2-4	w1
Timothy Grass	3-4	g6
Tomato	3-4	f25
Wheat	2-4	f4

Hitachi Chemical Diagnostics, Inc.
630 Clyde Court
Mountain View, California 94043
Tel. (650) 961-5501
Fax (650) 969-2745

Hitachi Chemical Diagnostics, Inc.
Hitachi Europe Ltd.
Whitebrook Park
Lower Cookham Road
Maidenhead, Berkshire SL6 8YA
United Kingdom
Tel. +44 (0) 1628 585 590
Fax. +44 (0) 1628 585 594

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