# HITACHI

 Hitachi Chemical Diagnostics, Inc.

## CLA<sup>®</sup> IgE ALLERGY NEGATIVE CONTROL REAGENT

Doc. No. 0441 - ENG  
Rev.: 12  
Lot No.: 1273-79  
Control Kit Part Number: 82056  
Bottle Part Number: 93523  
Control Kit Contents: 3 x 3 mL

### 1 Intended Use

CLA IgE Negative Control Reagent is intended for use in quality control procedures, which evaluate performance of the CLA Allergen-Specific IgE Assay.

### 2 Summary and Principle

The use of quality control materials is indicated as an objective assessment of assay performance.

### 3 Reagent

CLA IgE Negative Control Reagent is prepared from human serum. The control is provided in liquid form.

### 4 Storage and Stability

When stored at  $-20 \pm 10^{\circ}\text{C}$ , the CLA IgE Negative Control Reagent is stable until the expiration date indicated on the vial label. The reagent is supplied in volumes of 3 mL. When ready to use, thaw, remove the volume needed and re-freeze the unused material immediately. Do not store in the refrigerator. Unused portions of control reagents may be frozen at  $-20 \pm 10^{\circ}\text{C}$ . Only one freeze-thaw cycle is recommended. Thus aliquoting the control reagents into smaller volumes before refreezing may be helpful.

### 5 Expected Values

Each laboratory should establish its own mean values and acceptable ranges. Acceptable performance may be defined as a percentage of the results falling within the nondetectable range. Hitachi Chemical Diagnostics, Inc. defines acceptable performance criteria for this control as 85% of the allergens falling within the nondetectable class range.

### 6 Procedure

The CLA IgE Negative Control Reagent should be tested using the same procedure as patient specimens. Frequency of use may be determined by each Laboratory's Quality Control policies. Allow control reagent to equilibrate to room temperature. Invert gently several times and centrifuge for 10-20 minutes at 2000-3000xg or 2500-3600 rpm immediately prior to use.

### 7 Limitations

The CLA IgE Negative Control Reagent should not be used past the expiration date. If there is evidence of microbial contamination or excessive turbidity, discard the vial.

### WARNING


This product contains human source material. Treat as potentially infectious.

Each donor unit used to manufacture this control was tested by FDA accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Antigen, and antibody to HIV-1. No test method can offer complete assurance that products containing human source materials will be absent of these and other infectious agents. Handle this product with the same precautions used with patient specimens.

### FOR IN VITRO DIAGNOSTIC USE ONLY

Hitachi Chemical Diagnostics, Inc.  
630 Clyde Court  
Mountain View, California 94043  
Tel. (650) 961-5501  
Fax (650) 969-2745

Hitachi Chemical Diagnostics, Inc.  
Hitachi Europe Ltd.  
Whitebrook Park  
Lower Cookham Road  
Maidenhead, Berkshire SL6 8YA  
United Kingdom  
Tel. +44 (0) 1628 585 590  
Fax. +44 (0) 1628 585 594

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