OPTIGEN® POSITIVE CONTROL REAGENT

Doc. No. 0654 - ENG
 Rev.: 11
 Lot No.: 1244-120
 Control Kit Part Number: 82061
 Bottle Part Number: 93526
 Control Kit Contents: 3 x 3 mL

1 Intended Use

OPTIGEN Positive Control Reagent is intended for use in quality control procedures, which evaluate the performance of the OPTIGEN Assay.

2 Summary and Principle

The use of quality control materials is indicated as an objective assessment of assay performance. This positive control has been assayed for various specific allergens and their expected values are provided. Allergens that are not listed here may test positive with this lot of control reagent.

3 Reagent

OPTIGEN Positive Control Reagent is prepared from human serum with the addition of a preservative. The control is provided in liquid form.

4 Storage and Stability

When stored at $-20 \pm 10^{\circ}\text{C}$, the OPTIGEN Positive Control Reagent is stable until the expiration date indicated on the vial label. Once the reagent is defrosted, store the reagent at $2-8^{\circ}\text{C}$ for no longer than 8 weeks. Avoid freeze-thaw cycles.

5 Expected Values

The expected ranges printed in this insert were established by Hitachi Chemical Diagnostics, Inc. and are specific to this lot of OPTIGEN Positive Control Reagent. Individual laboratory means should fall within the expected range for each allergen. Procedural or reagent modifications may alter the mean values. Each laboratory should establish its own mean values and acceptable ranges and use those provided only as guides. Hitachi Chemical Diagnostics, Inc. defines acceptable performance criteria for this control as 85% of the assayed target allergens falling within the expected class range.

6 Procedure

OPTIGEN Positive Control Reagent should be tested using the same procedure as patient specimens. Frequency of use may be determined by each Laboratory's Quality Control policies. Allow control reagent to equilibrate to room temperature. Invert gently several times and centrifuge for 10-20 minutes at 2000-3000xg or 2500-3600 rpm immediately prior to use.

7 Limitations

The OPTIGEN Positive Control Reagent should not be used past the expiration date. If there is evidence of microbial contamination or excessive turbidity, discard the vial.

WARNING

This product contains human source material. Treat as potentially infectious.


Each donor unit used to manufacture this control was tested by FDA accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Antigen, and antibody to HIV-1. No test method can offer complete assurance that products containing human source materials will be absent of these and other infectious agents. Handle this product with the same precautions used with patient specimens.

FOR IN VITRO DIAGNOSTIC USE ONLY

Allergens	Expected Values Class Range	CLSI Code
American Sycamore	0	t11
Ash	1-3	t15
Aspergillus	2-4	m3
Baker's Yeast	2-4	f45
Beef	0	f27
Black Alder	0	t2
Casein	3-4	f78
Cat	2-4	e1
Cheddar Cheese	2-4	f81
Codfish	0	f3
Crab	1-3	f23
Cultivated Rye	2-4	g12
Egg White	2-4	f1
Hazelnut	1-3	f17
Housedust	2-4	h2
Latex	2-4	k82
Milk	2-4	f2
Mite Farinae	3-4	d2
Mite Pteronyssinus	3-4	d1
Peach	1-3	f95
Peanut	2-4	f13
Penicillium	2-4	m1
Salmon	0	f41
Shrimp	1-3	f24
Soybean	3-4	f14
Timothy Grass	2-4	g6
Wheat	2-4	f4

Hitachi Chemical Diagnostics, Inc.
 630 Clyde Court
 Mountain View, California 94043
 Tel. (650) 961-5501
 Fax (650) 969-2745

Hitachi Chemical Diagnostics, Inc.
 Hitachi Europe Ltd.
 Whitebrook Park
 Lower Cookham Road
 Maidenhead, Berkshire SL6 8YA
 United Kingdom
 Tel. +44 (0) 1628 585 590
 Fax. +44 (0) 1628 585 594

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