HITACHI

 **Hitachi Chemical Diagnostics, Inc.**

OPTIGEN® NEGATIVE CONTROL REAGENT

Doc. No. 0655 - ENG
Rev.: 05
Lot No.: 1273-21
Control Kit Part Number: 82062
Bottle Part Number: 93527
Control Kit Contents: 3 x 3 mL

1 *Intended Use*

OPTIGEN Negative Control Reagent is intended for use in quality control procedures, which evaluate performance of the OPTIGEN Assay.

2 *Summary and Principle*

The use of quality control materials is indicated as an objective assessment of assay performance.

3 *Reagent*

OPTIGEN Negative Control Reagent is prepared from human serum. The control is provided in liquid form.

4 *Storage and Stability*

When stored at $-20 \pm 10^{\circ}\text{C}$, the OPTIGEN Negative Control Reagent is stable until the expiration date indicated on the vial label. Once the reagent is defrosted, store the reagent at $2-8^{\circ}\text{C}$ for no longer than 8 weeks. Avoid freeze-thaw cycles.

5 *Expected Values*

Each laboratory should establish its own mean values and acceptable ranges and use those provided only as guides. Hitachi Chemical Diagnostics, Inc. defines acceptable performance criteria for this control as 85% of the allergens falling within the nondetectable class range.

6 *Procedure*

OPTIGEN Negative Control Reagent should be tested using the same procedure as patient specimens. Frequency of use may be determined by each Laboratory's Quality Control policies. Allow control reagent to equilibrate to room temperature. Invert gently several times and centrifuge for 10-20 minutes at 2000-3000xg or 2500-3600 rpm immediately prior to use. The control reagent should be run using the test procedure indicated for patient specimens.

7 *Limitations*

The OPTIGEN Negative Control Reagent should not be used past the expiration date. If there is evidence of microbial contamination or excessive turbidity, discard the vial.

WARNING


This product contains human source material. Treat as potentially infectious.

Each donor unit used to manufacture this control was tested by FDA accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Antigen, and antibody to HIV-1. No test method can offer complete assurance that products containing human source materials will be absent of these and other infectious agents. Handle this product with the same precautions used with patient specimens.

FOR IN VITRO DIAGNOSTIC USE ONLY

Hitachi Chemical Diagnostics, Inc.
630 Clyde Court
Mountain View, California 94043
Tel. (650) 961-5501
Fax (650) 969-2745

Hitachi Chemical Diagnostics, Inc.
Hitachi Europe Ltd.
Whitebrook Park
Lower Cookham Road
Maidenhead, Berkshire SL6 8YA
United Kingdom
Tel. +44 (0) 1628 585 590
Fax. +44 (0) 1628 585 594

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