

# Test Record

## CLA® Allergen-Specific IgE Assay

Patient Name

Patient ID

Test Date

Ordering Physician

Test Facility

Technician

LU	CLA® Class	Interpretation
> 242	4	Very High Positive
143 - 242	3	High Positive
66 - 142	2	Moderate Positive
27 - 65	1	Low Positive
12 - 26	1/0	Equivocal
0 - 11	0	No Allergy Detected

Test manufactured by:

# HITACHI

Hitachi Chemical Diagnostics, Inc.  
630 Clyde Court  
Mountain View, CA 94043-2239  
Tel: 800 233 6278  
Fax: 650 969 2745

[www.hcdiagnostics.com](http://www.hcdiagnostics.com)  
Doc. No. 0508 Rev. 06

16-JUNE-2008 14:25:18

DEVICE DIAGNOSTICS: PASS

\*-----\*

PATIENT ID:

PANEL: 84507  
Southeastern Inhalant

INTERNAL TEST CHECKS:  
Procedural >243 PASS  
Blanking 33 PASS

No. LU Allergen Class

1	90	Acacia	2
2	>300	Ash, White	4
3	49	Bch, Americn	1
4	39	Brch/Aldr Mx	1
5	72	BxEldr, Mple	2
6	77	Cedar, Mtn	2
7	2	Cottnwd, East	0
8	>300	Elm, White	4
9	133	Melaleuca	2
10	>300	Mulberry Mix	4
11	>300	Oak, White	4
12	>300	Pine Mix	4
13	108	Privet	2
14	10	Sycmr, Amrcn	0
15	0	Wlnt/Hck/Pcn	0
16	4	BahiaGrass	0
17	0	BermudaGrass	0
18	0	TimothyGrass	0
19	0	Cocklebur	0
20	4	Eng Plantain	0
21	2	Lamb's Qrtrs	0
22	1	Mrshldr, Rgh	0
23	0	Pigweed	0
24	5	Ragwd, Short	0
25	6	Sheep Sorrel	0
26	0	Waterhemp	0
27	2	Cat	0
28	3	Dog	0
29	0	CockroachMix	0
30	1	Mite, farinae	0
31	8	Housedust	0
32	0	Alternaria	0
33	4	Aspergillus	0
34	75	Candida	2
35	4	Cladosporium	0
36	0	Penicillium	0

END TEST

16-JUNE-2008 14:25:33